



epiCam® V  
Video Ophthalmoscope

veterinary technology for life  
**Eickemeyer®**

# All eyes on you

**NEW**



## epiCam® V Video Ophthalmoscope

epiCam® V is a high resolution video ophthalmoscope for veterinary use. View and capture live retinal video, then extract the best still images during review and document the images you require. The required software is included. Used with a standard Windows laptop, epiCam® V is lightweight, portable and non-invasive. Dilation is not mandatory, although it makes it slightly easier to see to the very wide periphery.


### Benefits

- ▶ Real-time, lossless retinal video examination
- ▶ Wide 45° field of view
- ▶ Premium optical quality and resolution
- ▶ Still images can be taken at up to 5 megapixels, ideal for publishing
- ▶ USB powered
- ▶ Patented design eliminates corneal glare

### Technical Specifications

- ▶ Interface: USB 2.0 or USB 3.0
- ▶ Operating systems supported: Windows 7, 8, 8.1, & 10
- ▶ Laptop processor required: Intel Core i5 for USB 3 mode
- ▶ Pixels: 1.296 x 972 pixels (USB 2.0); 2.592 x 1.944 pixels (USB 3.0)
- ▶ Field: 33° vertical x 45° horizontal
- ▶ Diameter of pupil: 4 mm or more
- ▶ Working distance: 13 mm
- ▶ Focal adjustment range: in excess of -15 D to +15 D
- ▶ Length: 153 mm
- ▶ Weight: 185 g
- ▶ Incl. footswitch

**173552**



### EICKEMEYER®

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
F 020 8891 2686  
E [info@eickemeyer.co.uk](mailto:info@eickemeyer.co.uk)

TELEPHONE 020 8891 2007

[www.eickemeyer.co.uk](http://www.eickemeyer.co.uk)