




PinCube

# *Pin them down!*

**NEW**



## **PinCube**

As part of an osteosynthesis procedure, bone and bone fragments are fixed for mutual stabilisation. This is usually done with plates, screws, tension bands and/or rods.

For some applications, such as tension band wiring, it is recommended that k-wires are positioned parallel to each other. This is free hand hardly to be managed.

For this reason, EICKEMEYER® has developed the PinCube!

### Advantages

- ▶ K-wires can be positioned exactly parallel
- ▶ The unique spiral-shaped hole pattern allows pins to be set both different distances and different diameters
- ▶ Usage of pins, drills and thread cutters in diameters of 0.6–2.0 mm
- ▶ Due to its extremely compact dimensions of maximum 35 mm, the PinCube is a very practical tool for small animal surgery

**191348**

#### **EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
F 020 8891 2686  
E [info@eickemeyer.co.uk](mailto:info@eickemeyer.co.uk)

**TELEPHONE 020 8891 2007**

**[www.eickemeyer.co.uk](http://www.eickemeyer.co.uk)**