



## Loupes & Loupe Light System

veterinary technology for life  
**Eickemeyer®**

# Take a Look!

## Loupes

The inexpensive magnifying glasses for demanding detailed work, convenient working distance and large visual field.

- ▶ 2x magnification
- ▶ Dioptre adjustment  $\pm 3$  dpt, right/left can be individually adjusted
- ▶ Excellent image quality through large focal area
- ▶ Supplied in a hard foam case and includes a lanyard for the glasses
- ▶ Weight: 40 g


**310600**

## Loupe Light System

The small, lightweight LED head lamp combines two high-power LEDs with a light condensing lens that produces a strong, brilliant and uniform light field.

- ▶ Allows a shadow-free and almost constant illumination for up to 8 hours in normal mode
- ▶ Vertically adjustable light
- ▶ Light Battery case (60 g) can be comfortably worn on the belt or in a shirt pocket
- ▶ Battery powered (3x AA batteries)

**310601**



### EICKEMEYER®

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
F 020 8891 2686  
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

[www.eickemeyer.co.uk](http://www.eickemeyer.co.uk)