



JUMPER Behandlungsstuhl
and 360° support for Jumper

veterinary technology for life
Eickemeyer®

I see your true colours



610155



JUMPER Saddle Stool

- ▶ Easily washable synthetic leather
- ▶ Available in various colours
- ▶ Smooth running double castors
- ▶ Seat width: 43 cm
- ▶ Adjustable height: 63–85 cm
- ▶ Aluminium base diameter: 54 cm

610155	Black	61015504	Anthracite Grey
61015501	Flame Red	61015505	Orange
61015502	Light Green	61015506	Pink
61015503	Night Blue		



610157

360° support for JUMPER

- ▶ 360° rotatable arm/backrest that supports one or both arms and can also be used as a lumbar support
- ▶ Height adjustable
- ▶ Optimum comfort with padded and seamless upper sides

610157

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
E info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk