



EICKEMEYER® LED  
Wood's Lamp

veterinary technology for life  
**Eickemeyer®**

# *Take a closer look*



## ***EICKEMEYER® LED Wood's Lamp***

The Wood's Lamp is a special diagnostic tool designed to illuminate and magnify skin, enabling clinicians to accurately inspect and diagnose various conditions. This lightweight and user-friendly LED lamp helps identify fluorescent diseased areas including bacterial infections, and pigment disorders.

- ▶ Window diameter: 65 mm
- ▶ 1.5x magnification
- ▶ Working distance: approx. 10 cm
- ▶ High depth of field
- ▶ 16 LEDs for uniform illumination
- ▶ Wavelength: 365 nm
- ▶ Weight: 305 g

**179030**

**EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**