



MAGIC 3000 PLUS
Portable Ultrasound Unit

veterinary technology for life
Eickemeyer®

The new lightweight for all who like to start with ultrasound

MAGIC 3000 PLUS Portable Ultrasound Unit

The MAGIC 3000 PLUS impresses with many new functions and features. Advanced technologies including Power Doppler (PDI) and PW Doppler help evaluate blood flow in different clinical conditions and the improved image processing provides outstanding image quality.



- ▶ Modern image optimization tools provide a very good picture
- ▶ Images and reports can be transferred directly to the computer or organised in the internal memory
- ▶ With iStation™, patient information is always at hand
- ▶ Custom buttons allow personalization
- ▶ LCD monitor: over 30 cm diagonal, tiltable
- ▶ Cine mode with memory function (5 sec.)
- ▶ Image export in Windows formats
- ▶ Image caption from library
- ▶ Measurement capabilities: distance, area, volume, circumference, angle, vet pregnancy calculation
- ▶ Quick start-up in approx. 17 seconds
- ▶ Long battery life of up to 2 hours (optional)
- ▶ Picture modes: B, B + B, M, B + M
- ▶ PDI (Power Doppler) and PW Doppler
- ▶ 8-fold TGC: balancing the signal strength of different layers
- ▶ iTouch™: automated gain setting
- ▶ iClear™: improves image resolution
- ▶ iZoom™: image representation filling the monitor
- ▶ iBeam™: improves image resolution in real time
- ▶ Internal SSD storage
- ▶ Two USB ports
- ▶ Outputs: Composite video, VGA, S-video
- ▶ Dimensions (in cm): W 30 x H 36 x D 16
- ▶ Lightweight: Only 5.5 kg

304045

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk



MAGIC 3000 PLUS
Portable Ultrasound Unit

veterinary technology for life
Eickemeyer®



303710



304021



303711



304022



304042

Optional Accessories

6.5 MHz Micro-Convex Probe

- ▶ 5.0 – 8.5 MHz

303710

7.5 MHz Linear T-Shaped Probe

- ▶ 5.0 – 10.0 MHz

303711

5.0 MHz Linear Rectal Probe

- ▶ 4.0 – 7.5 MHz

304021

7.5 MHz Linear Rectal Probe

- ▶ 5.0 – 10.0 MHz

304022

Battery

- ▶ Operation time: up to 2 hours

304041

Trolley

- ▶ Dimensions (in cm): W 47 x H 87 x D 54

304042

Carry Case

- ▶ Suitable for protecting the ultrasound unit

304043

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk