




Low Lift Exam/Prep Table

veterinary technology for life  
**Eickemeyer®**

# Goin' Down



Tie Down Cleat



## *Low Lift Exam / Prep Table*

The hydraulic lift mechanism helps protect your back. The heavy duty frame ensures stability throughout its range. All 4 castors can be locked simultaneously with the foot switch.

- ▶ Minimum height: 49,5 cm
- ▶ Maximum height: 98 cm
- ▶ Maximum load: 100 kg
- ▶ Four lockable castors
- ▶ Table Top dimensions (in cm): L 122 x W 60
- ▶ Includes: 6 tie-down cleats

**600200**

**EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
F 020 8891 2686  
E [info@eickemeyer.co.uk](mailto:info@eickemeyer.co.uk)

**TELEPHONE 020 8891 2007**

**[www.eickemeyer.co.uk](http://www.eickemeyer.co.uk)**