



Rabbit/Rodent Table Top
Mouth Gag

veterinary technology for life
Eickemeyer®

Safe Dental Treatments



Rabbit/Rodent Table Top Mouth Gag

The anaesthetised pet can be optimally fixed in order to gain the best possible access to the oral cavity. The worktable height and tilt are individually adjustable to achieve the optimal working position. The mouth opener is adjustable to fit any sized rabbit/rodent.

- ▶ Length (tabletop): 28 cm
- ▶ Width (between the struts): 13.5 cm
- ▶ Height: 29 cm

177324

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk