




RespiVet NANO
Respiratory Monitor

veterinary technology for life
Eickemeyer®

Every breath you take



RespiVet NANO Respiratory Monitor

Ultracompact and easy to use respiratory monitor that senses respiration from a 500 g kitten to a 100 kg Great Dane. Integrated alarm alerts you if the patient stops breathing for 30 seconds.

Technical Specifications

- ▶ Distinctive apnoea alarm
- ▶ Conveniently mounted on elbow adaptor
- ▶ On/Off switch
- ▶ Replaceable battery
- ▶ Patient Hose Connector: Outer Ø: 22 mm, Inner Ø: 15 mm
- ▶ Tracheal Tube Connector: Outer Ø: 15 mm

217110

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk