


SHOCKWAVE SYSTEM LIKAWAVE VARIO VET

Extracorporeal Pain Treatment


**Book your
demo today**



Equine

SHOCKWAVE SYSTEM LIKAWAVE VARIO VET

LIKAWAVE VARIO VET is a genuine innovation in high-energy shockwave therapy with its VARIO LOGIC TECHNOLOGY.



Effective

Vario Logic Technology intelligently modulates the shockwave frequency to allow cells to regenerate over and over again between impulses.

Gentle

Shockwave frequency and intensity (25 levels) can easily be adjusted to treat a wide range of conditions.

Targeted

Specific applicators deliver energy to the area depending on the indication, giving the desired depth of penetration.

Powerful

Maximum applied pressure of 0.58 mj/mm²

Extracorporeal Shockwave Therapy has been effective in the treatment of sport horses for many years.

User-friendly

Book your demo
T: 0208 891 2007 or
E: info@eickemeyer.co.uk

Eickemeyer®



SHOCKWAVE SYSTEM LIKAWAVE VARIO VET

- ▶ Unique VARIO LOGIC TECHNOLOGY
- ▶ 10.4" Touchscreen
- ▶ Pre-programmed VARIO Modes or own settings possible
- ▶ Various applicators, individually adjusted to the specific indications
- ▶ Quick access to easy-to-manage patient database
- ▶ 3 Different treatment modes: Linear, Triangular and Pulse Shape
- ▶ Optional Practice Trolley or Mobile Case/Trolley



Includes:

- | | | |
|--------------------------------------|--|-----------------------------|
| ▶ Handpiece with 3 m Connection Hose | ▶ Applicator S:
10 mm diameter – Green | ▶ Touchscreen Pen and Cloth |
| ▶ LIKAWAVE VARIO VET Accessory Case | ▶ Applicator M:
15 mm diameter – Blue | ▶ Ultrasound Gel (250 ml) |
| | ▶ Applicator L:
20 mm diameter – Purple | ▶ 3 m Power Cable |

326000

Equine

ACCESSORIES FOR LIKAWAVE VARIO VET

The following accessories are available to be purchased separately.



LIKAWAVE VARIO VET PRACTICE TROLLEY

- ▶ 4 large castors
 - ▶ Handpiece holder
 - ▶ Integrated Accessory Case
- 326003**



LIKAWAVE VARIO VET CASE

- ▶ With foam insert
 - ▶ Colour: Black
- 326001**



LIKAWAVE VARIO VET CASE TROLLEY

- ▶ Excludes Case
- 326002**

