




X-Ray Reference Ball

Always in the picture




X-Ray Reference Ball


The X-Ray Reference Ball on a flexible holder enables a scaled reference for X-Ray images, especially when assessing externally provided X-Ray images.

- ▶ Ø 25 mm
- ▶ Made of stainless steel
- ▶ Stable base with flexible neck
- ▶ Suitable for digital and analog systems

191990



application*



application 2**

*© Dr. Christoph Werner, Tierärztliche Gemeinschaftspraxis, Lohenstraße 5, 83395 Freilassing, Germany

**© Dr. Daniel Köch, Kleintierchirurgie AG, Ziegeleistraße 5, 8253 Diessenhofen, Switzerland