



EICKEMEYER® PinCube

veterinary technology for life
Eickemeyer®

Precise Pin Positioning

EICKEMEYER® PinCube

As part of an osteosynthesis procedure, bone and bone fragments are fixed for mutual stabilisation. This is usually done with plates, screws, tension bands and/or rods. For some applications, such as tension band wiring, it is recommended that k-wires are positioned parallel to each other. This is challenging to accurately do by hand.

For this reason EICKEMEYER® has developed the PinCube!

Advantages

- ▶ Enables K-wires can be positioned exactly parallel
- ▶ Unique spiral-shaped hole pattern allows pins to be set both different distances and different diameters
- ▶ Usage of pins, drills and thread cutters in diameters of 0.6 – 2.0 mm
- ▶ Extremely compact dimensions: 35 mm
- ▶ Very practical tool for small animal surgery

191348

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk