



RespiVet Respiratory Monitor

Gentle diagnostics



RespiVet Respiratory Monitor

The RespiVet Respiratory Monitor is an easy to use and low-cost monitoring solution. The unit has variable sensitivity settings and adjustable audible and visual apnoea alarms.

Technical Specifications

- ▶ Apnoea Alarm settings: 10–60 seconds
- ▶ Adjustable sensitivity for patients: 0,5 to 1.000 kg
- ▶ Photo eye for handsfree alarm reset
- ▶ LED readout of respiration rate
- ▶ Acoustic signal for each breath
- ▶ Connects endotracheal tube to the anaesthetic circuit
- ▶ Patient Hose Connector: Outer Ø: 22 mm, Inner Ø: 15 mm
- ▶ Tracheal Tube Connector: Outer Ø: 15 mm
- ▶ Rechargeable battery: Up to 70 hours of use
- ▶ Dimensions (in cm): L 13 x W 13 x H 6
- ▶ Includes: Connector (217106) and Sensor (217101)

217100

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk