



EickSorber
Anaesthetic Gas Adsorber

veterinary technology for life
Eickemeyer®

For your safety during anaesthesia



213488

EickSorber Anaesthetic Gas Adsorber

- ▶ Activated carbon filter to absorb organic anaesthetic gases
- ▶ Connects to pressure relief valve, thus preventing unwanted gases from entering room atmosphere.
- ▶ Replace filters after exceeding 1,000 g

Single Unit

- ▶ Without exhaust hose

21348801

Single Unit

- ▶ Includes 2 m exhaust hose (Ø 22 mm)

213488

Pack of 6 pieces

- ▶ Includes 2 m exhaust hose (Ø 22 mm)

21348806

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk