



EickTron 160 W
Electrosurgical Unit

veterinary technology for life
Eickemeyer®

Electric surgery made easy



EickTron 160 W ElectroSurgical Unit

The EickTron ElectroSurgery Unit is a versatile tool in veterinary medicine allowing it to be used for both monopolar and bipolar instruments. The unit's power adjustment is actuated via a foot or hand switch.

Includes

- ▶ Handpiece with finger switch and connecting cable
- ▶ Electrode set (10 pieces)
- ▶ Cable for neutral electrode
- ▶ Neutral electrode
- ▶ Footswitch
- ▶ Optional Bipolar accessories available

Specifications

- ▶ Input: 220–240 V; 50–60 Hz
- ▶ Power consumption: 300 VA
- ▶ Frequency: 600 kHz
- ▶ Max. Output 160 W (cut)
- ▶ Max. Output 100 W (coagulate)
- ▶ Max. Output 60 W (bipolar)
- ▶ Includes: Rectal probe
- ▶ Dimensions (in cm): H 10.4 x B 25.4 x T 28.8

323135

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk



EickTron 160 W Electrosurgical Unit

veterinary technology for life
Eickemeyer®



323105



323106



323139



323137



32313606



32313603




32313607



32313602



32313608



32313605



32313609



32313604



32313601



32313610



323138



323140

Accessories for EickTron Electrocautery Unit

Item No.	Description	Pieces
Bipolar Accessories		
323105	EickTron Bipolar Forceps • Autoclavable • Angled • Length: 150 mm	1
323106	EickTron Bipolar Forceps Cable • Autoclavable	1
323139	EickTron Bipolar Adaptor • Autoclavable	1
Monopolar Accessories		
323137	EickTron Monopolar Handpiece • Spare part	1
32313606	EickTron Electrode Angled Hook (L5)	5
32313603	EickTron Electrode (L2)	5
32313607	EickTron Electrode Needle - Angled (L6)	5
32313602	EickTron Electrode Ball - Angled (L10) • Ø 3.0 mm	5
32313608	EickTron Electrode Drop-shaped (L7)	5
32313605	EickTron Electrode Wire Loop (L4) • Ø 8.0 mm	5
32313609	EickTron Electrode Noose (L8)	5
32313604	EickTron Electrode Wire Loop (L3) • Ø 4.0 mm	5
32313601	EickTron Electrode (L1)	5
32313610	EickTron Electrode Ball - Straight (L9) • Ø 3.0 mm	5
323138	EickTron Neutral Electrode Cable • For monopolar electrodes	1
323140	Handle with Finger Switch	1
323141	EickTron Neutral Electrode	1
323143	Monopolar Scissors • Length: 18 cm	1
323144	Connection Cable • For Monopolar Scissors (323143) • Length: 3 m	1

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk