



EICKEMEYER®  
Ultrasonic Doppler

veterinary technology for life  
**Eickemeyer®**

# Accurate & reliable



217210

## **EICKEMEYER® Ultrasonic Doppler**

Reliable and accurate non-invasive blood pressure monitor for small animals.

- ▶ Compact hand-held design
- ▶ High sensibility for a clear signal
- ▶ Adjustable sound volume and headphone connection
- ▶ Focused 8 MHz Probe with wide beam for better flow detection in the artery
- ▶ Probe adapted for positioning on the carpal joint or the underside of the tail
- ▶ Rechargeable batteries are not interchangeable
- ▶ Ultrasonic frequency: 8 MHz
- ▶ Operating time: 6 hours
- ▶ Dimensions (in cm): L 17.5 x W 8 x H 2.5
- ▶ Weight: 300 g

**217210**



217308



217309

### **Sphygmomanometer**

- ▶ Easy one-hand operation
- ▶ Includes two cuffs: 5 cm and 7.5 cm

**217309**

### **Blood Pressure Cuff 3.3 cm**

- ▶ For cats and small dogs

**217308**

### **Blood Pressure Cuff 10 cm**

- ▶ For large dogs

**217311**

### **Blood Pressure Cuff 1.9 cm**

- ▶ For small cats

**217307**

**EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
F 020 8891 2686  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**