



SuperCut Scissors (MC)

veterinary technology for life
Eickemeyer®

The Razor's Edge



118311




118309



118310



118315



118315



118316

SuperCut Scissors (MC)

- ▶ Designed with micro-fine serrations on one blade to prevent soft tissue slipping while the opposing razor sharp blade cuts with minimal tissue trauma
- ▶ German stainless steel construction and blade design allows for a much longer working life before sharpening is required
- ▶ Ergonomically designed ring handles improve instrument control, provide a better feel when dissecting tissue and help reduce hand fatigue

- ▶ Straight
- ▶ Length: 90 mm

118309

- ▶ Curved
- ▶ Length: 90 mm

118310

- ▶ Straight
- ▶ Length: 115 mm

118311

- ▶ Curved
- ▶ Length: 115 mm

118312

SuperCut Metzenbaum Scissors (MC)

- ▶ Straight
- ▶ Length: 145 mm

118315

- ▶ Curved
- ▶ Length: 145 mm

118316

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk