



AlgerBrush II  
Corneal Burr

veterinary technology for life  
**Eickemeyer**®

# For Your Eyes Only

## AlgerBrush II Corneal Burr

The AlgerBrush II is a battery powered, handheld and low torque motorised burr. It has been designed to remove foreign bodies from the cornea or sclera of patients.

The burr in this rotary instrument is used as a “brush” rather than a drill and the foreign bodies are thus “brushed” from the eye leaving a smooth surface, which in turn heals much faster than with other instruments that could be used for this purpose.

According to the paper “Veterinary Ophthalmology” [(2013) 16, 2, 83-88], it has been successfully used to treat spontaneous corneal epithelial defects (SCCEDs) in dogs in combination with a bandage contact lens.

- ▶ Includes 3.5 mm round, medium grit diamond burr

**170010**

### AlgerBrush II Chuck and Burr

- ▶ 3.5 mm round, medium grit diamond burr

**170011**

### AlgerBrush II Corneal Burrs

- ▶ 1 mm
- ▶ 5 per pack

**170012**

### AlgerBrush II Chuck and Burr


- ▶ 2.5 mm round, medium grit diamond burr

**170013**

### AlgerBrush II Diamond Grinding Wheel (4 mm)

- ▶ Fine
- ▶ Convenient for recurrent corneal erosion

**170014**



**EICKEMEYER**®

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**