




EICKEMEYER® LED 4
Examination Light

veterinary technology for life
Eickemeyer®

Let there be light!



EICKEMEYER® LED 4 Examination Light

Modern examination light for the practice with energy-saving LED technology. With an illuminance of approx. 50,000 lux at a working distance of 30 cm, the EICKEMEYER® LED 4 examination light offers optimal illumination for a professional examination.

Technical specifications

- ▶ 15 V / 11.25 W LED lamp with 4 high-performance LEDs
- ▶ Illuminance: approx. 50,000 lux / 30 cm distance
- ▶ Aluminum light head for optimal heat convection and heat radiation
- ▶ Lifespan: approx. 50,000 hours
- ▶ Color temperature: approx. 4,200 K
- ▶ Gooseneck construction with 90 cm long arm (arm remains in the desired position)
- ▶ Maximum examination height: 165 cm
- ▶ Stable, tilt-proof mobile stand (5-foot stand) made of powder-coated metal

611111

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk