



EickWrap Cohesive Bandages

veterinary technology for life
Eickemeyer®

Strength and lightness combined



050750

050751

050752



EickWrap Cohesive Bandages

- ▶ Highest quality cohesive bandage
- ▶ Conforms itself to the shape of the area where it is applied – helping to keep dressing in place
- ▶ Self-adhesive but, it will not stick to the skin
- ▶ Contains latex
- ▶ Allows air to circulate over the skin and around the dressing
- ▶ Individually packaged – assorted colours
- ▶ Length: 4.5 m / roll

- ▶ Width: 5 cm
- ▶ 36 per pack

050750

- ▶ Width: 7.5 cm
- ▶ 24 per pack

050751

- ▶ Width: 10 cm
- ▶ 18 per pack

050752

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk