



HiRay 5 X-Ray Machine

Straight to the bone



701105

HiRay 5 X-Ray Machine

The portable HiRay 5 combines ultimate power with extreme compactness. It produces a nominal output of 5 kW, while weighing only 14 kg! The greater mA range allows for shorter exposure times. This is ideal for imaging regions with rapid motion such as the small animal thorax.

- ▶ High-voltage setting: 120 kV/100 mAs (30 mA)
- ▶ High-current setting: 50 kV/32 mAs (100 mA)
- ▶ kV setting: 40 to 120 kV (in 1 kV steps)
- ▶ mAs settings: 0.3 to 100 mAs (26 settings)
- ▶ Generator: High frequency inverter, 5 kW nominal output
- ▶ X-Ray Tube Focal spot: 1.8 mm
- ▶ Collimator: Field lamp with a 30 second timer and across indication
- ▶ Line cable: 3 metres
- ▶ Hand switch: Two stage, deadman type with curled cable
- ▶ Dimensions (in cm): L 35 x W 20.3 x H 16
- ▶ Weight: 14 kg
- ▶ Foot switch optional

701105

Portable X-Ray Stand

701118

HiRay Foot Switch

701121



701118



701121

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk