




Disposable Surgical Drape
and Wall Dispenser for
Surgical Drape

veterinary technology for life
Eickemeyer®

Clean as a whistle




920060

Disposable Surgical Drape

- ▶ Ideal for covering the whole operating table
- ▶ Simple way to reduce germs in the operating theatre
- ▶ Double folded
- ▶ Thickness: 0.025 mm
- ▶ Dimensions: W 120 cm x L 200 m

920060



920062

Wall Dispenser for Surgical Drape

- ▶ Simple, single handed operation
- ▶ Spring loaded tear-off rail for constant pressure and good separation
- ▶ Easy roll change with hinged tear-off rail
- ▶ For drape rolls with a width of 60 cm, Ø 24 cm and max. 20 kg
- ▶ Dimensions (in mm): W 666 x D 215 x H 97

920062

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk