



EickView Ultra Compact  
Endoscope Set

veterinary technology for life  
**Eickemeyer®**

# *The new generation of rigid endoscopes*

## ***EickView Ultra Compact Endoscope Set***

The EickView Ultra Compact is the new generation of multifunctional rigid endoscopes. It is ideal for all ENT applications as well as tracheo- and vaginoscopies with its 3.0 mm diameter.

### **Features**

- ▶ Homogeneous and minimally invasive distal end by fusion of workforce and endoscope
- ▶ No dirt deposit due to homogeneous distal end
- ▶ Slightly smooth access due to extremely low outer diameter
- ▶ Integrated working channel allows for the removal of foreign body extractions and biopsies

### **Technical specifications**

- ▶ Outside Ø: 3 mm
- ▶ Working length: 11 cm
- ▶ Working channel: Ø 1.3 mm
- ▶ Viewing angle: 0°
- ▶ Opening angle: 90° (wide-angle)
- ▶ Endoscope can be inserted
- ▶ Flushing taps: Two integrated flushing taps allow continuous insufflation and supply of flushing fluid
- ▶ Set includes: Endoscope, LED Light Source incl. charger, Spare battery, Foreign Body Forceps, 2 Cleaning brushes, Transport Case

**306112**

#### **EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
F 020 8891 2686  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**