



AESCULAP Isis Cat Clipper

veterinary technology for life
Eickemeyer®

Savor the silence



902605



902601



902513

AESCULAP Isis Cat Clipper

The AESCULAP Isis cordless clipper is ideal for preparation for small areas (forelimb, head and neck) when taking bloods, biopsy or for hard-to-reach areas.

- ▶ Quiet and low vibration
- ▶ Ergonomic design
- ▶ Cutting width: 24 mm
- ▶ Cutting length: 0.5 mm
- ▶ Li-ion battery: runtime up to 120 min with approx. 90 min charging time
- ▶ Charge indicator via charging station
- ▶ Weight: 120 g (incl. battery & blade)
- ▶ Includes: clipper, blade (carbon blade – cannot be sharpened), charging station, attachment comb (adjustable to 4 different cutting lengths, 3 – 6 mm), brush, AESCULAP oil, instructions

902605

Isis/Exacta Clipper Blade

- ▶ Carbon blades (can not be resharpened)

902601

EICKEMEYER® Clipper Spray

- ▶ Clipper blades require cleaning in-between each patient to provide the longest possible working life
- ▶ Dual action of disinfection and lubrication

902513

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk