



CT, MRI and Ultrasound
Positioner

veterinary technology for life

Eickemeyer®

Perfect position – perfect vision

CT, MRI and Ultrasound Positioner

- ▶ Facilitates positioning of patients to allow abdominal ultrasonography
- ▶ Perfect positioning helps keep the patient calm and comfortable
- ▶ Robust, easy to clean and disinfect

Size 1

- ▶ Suitable for cats and small dogs
- ▶ Dimensions (in cm): W 25 x H 10 x L 80

303991

Size 2


- ▶ Suitable for medium-sized dogs
- ▶ Dimensions (in cm): W 35 x H 14 x L 100

303992

Size 3

- ▶ Suitable for large-sized dogs
- ▶ Dimensions (in cm): W 50 x H 20 x L 130

303993



303992

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk