


Disposable Ring Retractor
with elastic hooks

veterinary technology for life
Eickemeyer®

A round thing



Disposable Ring Retractor with elastic hooks

A Ring Retractor facilitates access to the surgical field during various soft tissue surgery procedures.

- ▶ Four hinges allow adaptation to the body shape
- ▶ Makes an area of 14.4 x 14.4 cm
- ▶ Incl. 8 elastic hooks (sharp)
- ▶ Packaged sterile

162150

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk