



HAICO PRO 3
Equine Operating Table


veterinary technology for life
Eickemeyer®

The Ultimate Tool

HAICO PRO 3 Equine Operating Table

- ▶ Tilt and swivel adjustment in 4 directions
- ▶ Wider and longer load bearing surface
- ▶ Larger recess in the middle section allows more freedom of movement for the surgeon
- ▶ All supporting parts areas are easy to remove
- ▶ Head section can be moved laterally
- ▶ Maximum load 1,300 kg
- ▶ Dimensions (in cm): L 310 x W 145
- ▶ Weight: 500 kg

608760



EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk