



Portable Suction Pump

# Clear vision in the wound area



306050

## Portable Suction Pump

This unit is designed to combine both functionality and aesthetic appeal, while adhering to today's strict standards of hygiene and infection control. This product is intended for use in multiple areas of the practice – operating theatre, prep room or in consulting rooms.

- ▶ Includes one filter
- ▶ Suction 40 l/min
- ▶ Vacuum: -0.40 to -0.80 bar
- ▶ Container: 1 litre
- ▶ Dimensions (in cm): L 35 x W 21 x H 18
- ▶ Weight: 3.6 kg

**306050**

## Accessories

Suction Tip, curved,  
L 17 cm, Ø 5 mm

**199219**

FERGUSON Suction Tip,  
curved, L 19 cm, Ø 3 mm

**199230**

YANKAUER Suction Tip,  
curved, L 28.5 cm, Ø 2 mm

**199220**

POOLE Suction Tip, curved,  
L 22 cm, Ø 8 mm


**199225**

FERGUSON Suction Tip,  
curved, L 19 cm, Ø 5 mm

**199250**

YANKAUER Suction Tip,  
curved, L 28.5 cm, Ø 4.8 mm

**199221**



199219

199225

199250

199220

### EICKEMEYER®

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**