



OxyVet III
Oxygen Concentrator

Oxygen 4 U



213175

OxyVet III Oxygen Concentrator

An oxygen concentrator draws in room air, separates the oxygen from other gases in the environment and delivers the oxygen at high concentrations.

Features

- ▶ Flowrate: 0,5 to 5 l/min
- ▶ Oxygen concentration: 93 % ± 3 %
- ▶ Sound level: 40 dbA (50 Hz)
- ▶ Electrical requirements: 230 V/50 Hz
- ▶ Power consumption: approx 290 watt @ 2 l/min and approx. 312 watt at 5 l/min
- ▶ Includes 2,1 metre Oxygen Hose
- ▶ Dimensions (in cm): L 34 x W 30,4 x H 62
- ▶ Weight: 16,3 kg
- ▶ 3 Year Warranty

213175

Bubble Humidifier for Oxygen Concentrators

- ▶ Humidifies oxygen prior to patient delivery
- ▶ Compatible with OxyVet III (213175) and OxyVet III PLUS (213177)

213176

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk