



HiRay PLUS X-Ray Machine

veterinary technology for life
Eickemeyer®

Seeing is Believing



701115

HiRay PLUS X-Ray Machine

One of the smallest and most powerful portable High Frequency X-Ray units in the world! The HiRay PLUS is a portable unit, which is suitable for outdoor use, as well as inside small animal practices.

Technical Specifications

- ▶ High RF power: 40 to max. 100 kV/30 mA, 2.2 kW
- ▶ 32 mAs settings, 0.3 to 20 mAs with illuminated LCD display
- ▶ Anatomical programming of 8 kV/mAs combinations
- ▶ Focus: 1.2 x 1.2 mm
- ▶ Collimator with bright halogen bulb and automatic timer
- ▶ Integrated laser beam for precise distance measurement
- ▶ 2-Speed manual exposure switch on 2.5 m curled cable
- ▶ Includes: practical and sturdy aluminum case
- ▶ CE Certified & PTB approved and tested
- ▶ Dimensions (in cm): L 28 x W 16.3 x H 16
- ▶ Weight: 8,9 kg

701115

Portable X-Ray Stand

701118

HiRay Foot Switch

701121



701118



701121

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk