



Speculum WORLD WIDE  
EQUINE 2020 Series

veterinary technology for life  
**Eickemeyer®**

# *For a safe dental treatment*



## ***Speculum WORLD WIDE EQUINE 2020 Series***

This speculum can be used on horses and cattle. When using it on cattle place it upside down and use a gum bar on top for best results.

- ▶ Double stitched leather straps for a secure and durable grip
- ▶ Stainless steel construction
- ▶ 15 adjustable height levels
- ▶ Improved ratchet mechanism
- ▶ Interchangeable latch bars
- ▶ Very easy to open
- ▶ Large side access

**177529**

**EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**