



CryoIQ® Cryosurgical Unit

veterinary technology for life
Eickemeyer®

Cold as ice



173443



173446

CryoIQ® Cryosurgical Unit

Cryosurgery (cryotherapy) is the local application of extreme cold for the controlled destruction or removal of abnormal or diseased tissue. The compact, portable and hand-held CryoIQ® Cryosurgical Device uses N₂O as a cryogen with a freezing temperature of up to -89 °C.

- ▶ Pencil handle for easy and very precise application
- ▶ Longer cartridge life thanks to the fine metal cap
- ▶ Almost painless and bloodless procedure
- ▶ Ideal for ophthalmic or dermatological treatments
- ▶ The degree of cold is determined by how much the lever is pressed

CryoIQ® DERM plus Cryosurgical Unit

- ▶ CryoIQ® DERM plus enables microcapillaries to freeze the area to be treated with millimeter precision.
- ▶ Includes: attachment for spray freezing (Ø 1–6 mm), gas cartridge N₂O (25 g), gauge, storage box, instructions for use

173443

CryoIQ® PRO Cryosurgical Unit

- ▶ CryoIQ® PRO offers a precise treatment of affected skin areas without damaging the surrounding healthy tissue.
- ▶ Advantages: More suitable for larger lesions, interchangeable attachments, higher freezing effects, longer tips make it easier to treat animals with more hair
- ▶ Includes: Standard attachment for spray freezing (Ø 1–6 mm, L 10 mm), dermatology attachment for spray freezing (Ø 7–18 mm, L 22 mm), PRO housing with lever, gas cartridge N₂O (25 g), gauge, storage box, instructions for use

173446

N₂O Cartridge for CryoIQ®

- ▶ Contents: 16 g
- ▶ Contents: 25 g

173444

173445

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
F 020 8891 2686
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk