

## STANDARD / DTE Scaling Tips

Item no. 176871 – 176874, 176878 / 17402401 – 17402405

### Differentiation between Scaler Tips supplied with EICKEMEYER® Dental Scalers

EICKEMEYER® sells two different types of scaler tips that are **not compatible** with each other. Using the wrong scaler tip will damage the thread on the handpiece and the handpiece will have to be replaced. Please ensure that you always use the correct scaler tips.

#### STANDARD Scaler Tips


- ▶ For handpiece with **red ring**
- ▶ For Profident LUMOS built-in scalers (17402319 – models before 04/2021)
- ▶ For EickSonic II Dental Scaler (176865)
- ▶ Compatible with Woodpecker and EMS

**176871 – 176874, 176878**


#### DTE Scaler Tips

- ▶ For handpiece with a **black ring**
- ▶ For Profident LUMOS built-in scalers (17402324 – models after 04/2021)
- ▶ The tips are imprinted with DTE
- ▶ Compatible with Satelek and NSK

**17402401 – 17402405**



STANDARD Scaler Tips



DTE Scaler Tips