





Vacuum Pillows &
Vacuum Pump

veterinary technology for life
Eickemeyer®


Safe positioning



605945



605980



605990

Vacuum Pillows

Vacuum pillow designed to help position and support an animal during surgical procedures. Adjusting the inflation pressure allows the pillow to adapt to a wide range of patient sizes.

- ▶ Constructed from durable PVC with extra soft edges
- ▶ Helps protect animals from hypothermia

Dimensions (in cm):

L 60 x W 30

605935

Dimensions (in cm):

L 100 x W 50

605945

- ▶ Four handles to carry the animal along with three velcro straps to help secure the animal
- ▶ Dimensions (in cm): L 100 x W 80

605980

Vacuum Pump

- ▶ Suitable for inflating all EICKEMEYER® vacuum cushions
- ▶ Robust, impact-resistant plastic

605990

EICKEMEYER®

3 Windmill Business Village
Brooklands Close, Sunbury-on-Thames
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

www.eickemeyer.co.uk